

PREPARING FOR THE ARTAS® ROBOTIC PROCEDURE — TRANSITIONAL HAIRCUTS

We encourage a transitional haircut starting a month or so prior to your ARTAS® Robotic Procedure. Ultimately, it is best to have your hair 1 millimeter in length the day of your procedure. Typically within 7 to 10 days after your ARTAS Robotic Procedure, your donor area will be nearly undetectable.

Steps for a Transitional Haircut



Start with your original hair cut



A month prior to your ARTAS Robotic Procedure, shorten your hair length to 3 millimeters



2-3 days prior to your ARTAS Robotic Procedure, shorten your hair length to 1 millimeter